



A highly configurable platform for advanced PPG analysis






Flavie Durand de Gevigney, Jean-François Nezan, Eric Bazin, Julien Heulot and Mickael Dardaillon.








Detecting change in volume using lights







Transmitted Reflected LED LED PD









List of possible measures and application of the PPG

Measures

Heart Rhythm

Heart Rhythm Variation

Oxygenation

Breathing

Blood pressure

Glycemia

Hydration

Application

Training monitoring

Stress monitoring

Cardiovascular health assessment

Training overload detection

Disease detection


Epileptic seizure detection








IETR Physical principal of the calculation of the oxygenation






Absorbance Molar absorption coefficient M-1cm-1 Molar concentration M optical path length cm


 $A = \varepsilon c l$



Classical sport watch



Proposed PPG sensor








UIFTR Our proposed platform



Photodiode and **LEDs**

Battery

Reset Button



Indicator LEDs

Control Button







- 8 wavelengths
- Highly configurable
 - Sampling frequency
 - Illumination parameters
 - etc
- High accuracy
- Choice between reflectance and transmittance
- Access to raw data








IETR New possibility with our proposed platform

- Be able to test different and new scenarii
- Determine optimal parameters
- Test energy aware scenarii
- Use in different light environment
- Adapt to the physiological differences between subjects






















More detail on the Led power stage
















More detail on the Led power stage










More details on the acquisition stage

How?













More details on the acquisition stage

How?





IETR <u>Differences between Analog and Digital Zoom</u>






AC + DC





Finding breathing pattern: proof of concept











IETR Future





Creation of a dataset

- Test of our sensor against different type of other sensors
- Large panel of activities

New applications

- Glycemia estimation
- SpO2 estimation
 - Dyshemoglobin estimation





IIIETR Conclusion

PPG in general

- Non invasive
- Easy to transport
- Easy to use
- Lot of measures possibles
- Lot of applications

Addition of our platform

- 8 wavelengths
- Highly configurable
- High accuracy
- Choice between reflectance and transmittance
- Access to raw data
- New measures
- New applications



